Q-PLUS™-Green Real Time qPCR Master Mix (2x)

Ordering Info

TBK0055, 5 reactions (sample)

TBK0057, 400 reactions

TBK0059, 1000 reactions

TBK0056, 200 reactions

TBK0058, 500 reactions

Description

Q-PLUS™-Green Real Time qPCR Master Mix (2X) is an optimized formula dye-based Master Mix to be used in real time quantitative PCR. This 2x reaction mixture contains all the components necessary for real-time PCR, except for primers and the template, enabling the amplification and detection of DNA. This kit also includes a separate vial of ROX that can be optionally added to the qPCR reaction Mix. The final concentration of ROX will vary depending on each real-time cycler manufacturer's specification.

Features

- Ready to use.
- Enhanced specificity and sensitivity, amplifies low copy number targets with reduced non-specific.
- Compatible with fast and standard PCR program.

Kit Components

Components	TBK0056	TBK0057	TBK0058	TBK0059
Q-PLUS™-Green Real Time qPCR Master Mix	2x 1mL	4x 1mL	5x 1mL	10x 1mL
PCR Grade Water	2 mL	2x 2mL	3x 2 mL	5x 2mL
ROX reference dye	1 vial	1 vial	1 vial	2 vial

Order Info Kit Components: ROX Reference Dye (TBR0278) | PCR Grade Water, nuclease free (TBB0303).

Storage

Shipped on blue ice. Upon receipt, kit components should be immediately stored at -20°C. Avoid repeated freezing and thawing. Maintain cold when thawed.

Applications

- Real time quantitative PCR to quantify DNA/ cDNA: gene copy number determination, microbial detection, gene expression.
- Genotyping.

Technical Assistance

Please refer any technical questions to support@tiarisbiosciences.com

ROX Reference Dye

The passive reference dye ROX is necessary for certain real-time PCR machines as it compensates for non-PCR-related variations in fluorescence detection. The fluorescence emitted from the ROX dye remains constant throughout the real time PCR process, providing a stable baseline against which PCR-related fluorescent signals can be normalized. As a result, the ROX dye can compensate for any differences in fluorescence detection between wells that may arise from slight variations in reaction volume or differences in well position.

Depending on your equipment, prior to use for the first time, add 2 µL (Low ROX) or 20 µL (High ROX) of the ROX reference dye to the 1 mL Q-PLUS™-Green Real Time qPCR Master Mix, (2x) and vortex briefly. Once ROX has been added, the Master Mix can be used directly or stored at -20°C for up to 1 year.

PROTOCOL

1. Gently vortex and briefly centrifuge kit components after thawing.
2. Place a tube on ice and add the following components for each 20 μ L reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras:

Components	Volume	Final Concentration
Q-PLUS™ Green qPCR Master Mix (2x)*	10 μ L	1x
Forward primer (5 pmol/ μ L)	<2 μ L	50-400 nM
Reverse primer (5 pmol/ μ L)	<2 μ L	50-400 nM
Template DNA (step 5)	**	**
PCR Grade Water	up 20 μ L	
Final Volume	20 μL	

* ROX as required; see section: "ROX Reference Dye"

** Optimal amounts of template: genomic DNA <1 μ g, plasmid or viral DNA 1-50 ng, cDNA synthesis reaction 1-5 μ L (or <100 ng).

3. Dispense the master mix into wells of PCR plate.
4. Gently vortex and spin down the samples.
5. Add in each well the DNA sample. Mix well by pipetting.
6. Seal the PCR plate with optical film.
7. Set-up qPCR cycling (if applicable, select fast mode on the instrument*):

Suggested thermal cycling conditions:

Process	Cycles	Temperature	Time
Enzyme Activation	1 x	95 °C	2-3 min
Denaturation	40 x	95 °C	5 sec
Annealing/ Extension		60-65°C	20-30sec**
Dissociation/Melt Analysis	See your instrument guidelines for setup.		

* For efficient qPCR, under fast cycling conditions, it is recommended to amplify DNA fragments ranging from 80-200bp.

** Select the shortest time possible but not less than 20 sec and do not exceed 30 seconds. Primers should have melting temperature of approximately 60°C.

In case of 3-step cycling, anneal at optimal annealing temperature for 20 sec and minimum time necessary at 72°C for data acquisition (according to manufacturer's guidelines).

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