High-Q™ -Spin-Columns Blood & Cell Culture Genomic DNA Purification Kit

Ordering info

TBK0121, 3 reactions (sample)

TBK0123, 200 reactions

TBK0122, 50 reactions

Description

High-Q™-Spin-Columns Blood & Cell Culture Genomic DNA Purification Kit is a silica-membrane-based DNA purification kit. to obtain genomic DNA with high quality and purity. Suitable for blood, plasma, serum and other body fluids as well as for cell culture samples.

Features

- **High yield and purity**, 2-15 µg, $A_{260}/A_{280} = 1.8 \pm 0.2$; $A_{260}/A_{230} = 2.0 \pm 0.2$.
- Suitable for **wide variety of samples** like fresh or frozen anticoagulated mammalian blood (treated with EDTA, Citrate or Heparine), aged blood, plasma, serum.
- Complete Blood Sample from **25 µL to 200 µL**.

Applications

DNA obtained is suitable for downstream molecular biology applications such as real-time PCR, endpoint PCR, multiplex PCR and enzymatic digestion for cloning or Southern, genotyping, etc.

Quality Control

DNA isolation from 200 µL whole human blood is checked by: integrity (agarose gel electrophoresis), quantity and quality ($A_{260}/A_{280} = 1.8 \pm 0.2$).

Material required (not supplied)

- Ethanol (CAS 64-17-5).
- Isopropanol (CAS 67-63-0)

Kit Components

Components	TBK0122	TBK0123
High-Q™ Spin Column with Collection Tubes	50	200
Blood Lysis Buffer	20 mL	65 mL
PBS 1x pH 7.4	15 mL	60 mL
Blood Binding Buffer	8 mL ^a	26 mL ^b
Blood Washing Buffer	12 mL ^c	44 mL ^d
WB2 Buffer	16 mL ^e	2 x 30 mL ^f
Elution Buffer	15 mL	25 mL
Proteinase K	30 mg ^g	4 x 30 mg ^g
Proteinase K Resuspension Buffer	1.5 mL	4 x 1.5 mL

Order Info Kit Components: High-Q™ Spin Column with Collection Tubes (TBM0010) | Blood Lysis Buffer (TBB0498) | PBS (TBB0360) | Proteinase K (TBZ0305) | Proteinase K Resuspension Buffer (TBB0546) | Blood Binding Buffer (TBB0499) | Blood Washing Buffer (TBB0500) | WB2 Buffer (TBB0512) | Elution Buffer (TBB0510).

Before its use:

- ^a Add 12 mL isopropanol and mix well.
- ^b Add 39 mL isopropanol and mix well.
- ^c Add 18 mL isopropanol and mix well.
- ^d Add 66 mL isopropanol and mix well.
- ^e Add 64 mL absolute ethanol and mix well
- ^f Add 120 mL absolute ethanol and mix well.
- ^g Add 1.5 mL Proteinase K Resuspension Buffer and mix well. Store at -20°C.

Storage


Store the kit at 25°C.

Store Proteinase K at -20°C.

PROTOCOL FOR GENOMIC DNA PURIFICATION

I. FROM BLOOD

1. Transfer **25 μ L** Proteinase K to a 1.5 mL tube (not provided).
2. Add to up **200 μ L** Sample.
3. *Optional: add 10 μ L RNase (20 mg/mL) (not provided).*
4. Add **300 μ L** Blood Lysis Buffer and mix by vortex until homogenous solution is observed.
5. Incubate at **60°C**, 15 minutes and shake occasionally.
6. Add **300 μ L** Blood Binding Buffer and mix by vortex vigorously.
 - ✓ Check isopropanol has been added to Blood Binding Buffer.
7. Place a High-Q™ Spin Column into a Collection Tube and transfer lysate from Step 6 to the spin column reservoir using a pipette.
8. Centrifugate at **13,000 g** for 1 minute and discard the flow-through.
9. Add **500 μ L** Blood Washing Buffer.
 - ✓ Check isopropanol has been added to Blood Washing Buffer.
10. Centrifugate at **13,000 g** for 1 minute and discard the flow-through.
11. Add **700 μ L** WB2 Buffer.
 - ✓ Check absolute ethanol has been added to WB2 Buffer.
12. Centrifugate at **13,000 g** for 1 minute and discard the flow-through. Repeat steps 11 and 12 one more time.
13. Centrifugate again at **13,000 g** for 2 minutes to remove residual ethanol.
14. Place the reservoir into a clean 1.5 mL tube.
15. Add **50-100 μ L** EB Buffer or Water, Molecular Biology Grade (pre-heated at 70°C) on top of the silica membrane.
16. Incubate at room temperature for 2 minutes.
17. Centrifugate at **13,000 g** for 1 minute to collect DNA in eluate.

II. FROM CULTURED CELLS

1. Transfer **5x10⁶ cells** to a 1.5 mL tube.
2. Centrifuge at **3,000 g** for 5 minutes and discard supernatant.
3. Add **200 μ L** PBS 1x pH 7.4 to the cell pellet and resuspend the pellet.
4. In a tube containing **25 μ L** Proteinase K, add **200 μ L** resuspended cells from previous step. Continue at **Step 3** of the previous procedure.

High-Q™ -Spin-Columns Blood/ Cell Culture Genomic DNA Purification Kit

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Referencia Componentes: High-Q™ Spin Column with Collection Tubes (TBM0010) | Blood Lysis Buffer (TBB0498) | PBS 1x pH 7,4 (TBB0360) | Proteinase K (TBZ0305) | Proteinase K Resuspension Buffer (TBB0546) | Blood Binding Buffer (TBB0499) | Blood Washing Buffer (TBB0500) | WB2 Buffer (TBB0512) | Elution Buffer (TBB0510).

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