

High-Q™ Spin-Column Gel Extraction & Cleanup Kit

Ordering info

TBK0190, 3 reactions (sample)

TBK0191, 50 reactions

TBK0192, 200 reactions

Description

High-Q™ Spin-Column Gel Extraction & Cleanup Kit is an easy silica-membrane-based DNA purification kit suitable for the purification of DNA embedded in agarose gels or to remove salts, proteins, nucleotides, primers or other contaminants from enzymatic reactions.

Features

- Easy and fast protocol.
- Complete removal of DNA contaminants.
- Excellent DNA recovery.
- No phenol extraction.

Applications

- DNA isolation from DNA embedded in agarose gels or enzymatic reaction mixtures such as PCR, restriction digestion, labelling reactions, dephosphorylation, etc.
- DNA obtained is suitable for downstream molecular biology applications such as cloning, PCR, sequencing, digestion, genotyping, etc.

Quality Control

Isolation of a 0.5 kb DNA fragment from 0.7% agarose gel is checked by: integrity (agarose gel electrophoresis), quantity and quality ($A_{260}/A_{280} = 1.8 \pm 0.2$).

Kit Components

Components	TBK0191	TBK0192
High-Q™ Spin Columns	50	200
Collection Tubes	50	200
BQ Buffer	60 mL	240 mL
BCL1 Buffer	30 mL	120 mL
BQ2 Buffer	8 mL ^a	35 mL ^b
Sodium Acetate 3M, pH 5	1 mL	1 mL
Elution Buffer	15 mL	15 mL
COBAL™ DNA Loading Buffer 6x	0.1 mL	0.5 mL

Order Info Kit Components: High-Q™ Spin Columns (TBM0017) | Collection Tubes (TBM0020) | BQ Buffer (TBB0527) | BCL1 Buffer (TBB0535) | BQ2 Buffer (TBB0529) | Elution Buffer (TBB0510) | COBAL™ DNA Loading Buffer 6x (TBB0321).

Before its use:

- ^a Add 32 mL absolute ethanol and mix well.
- ^b Add 140 mL absolute ethanol and mix well.

Storage

Store the kit at 25°C.

Material required (not supplied)

- Ethanol (CAS 64-17-5).
- Isopropanol (CAS 67-63-0).
- 1.5 mL Tubes.

PROTOCOL

I. DNA PURIFICATION FROM AGAROSE GELS

1. Cut with a clean scalpel the agarose region containing your DNA fragment. Try to eliminate all the extra agarose gel.
2. Introduce the slice gel into a pre-weighed 1.5 mL tube. Weigh the gel excised.

If the weight of gel excised is higher than 400 mg, split it in two pieces and use two columns.

3. Add BQ Buffer for each 100 mg of gel weight as it is indicated:

Agarose Gel Concentration	BQ Buffer Volume	Examples
< 2%	300 µL	To 110 mg agarose gel, add 330 µL BQ Buffer
≥ 2- 3%	500 µL	To 110 mg agarose gel, add 550 µL BQ Buffer

4. Incubate at 55-60 °C and shake occasionally until agarose is completely dissolved (aprox. 10 minutes).
5. *Optional, if your DNA fragment is <500 bp or >10 kb, add 1 gel volume of isopropanol and mix well.*
6. To guarantee optimal DNA adsorption, add 3-5 µL 3M Sodium Acetate pH=5.0 and mix well.
7. Transfer the mix to a High-Q™ Spin Column placed into a Collection Tube.

The maximum loading volume of High-Q™ Spin Column is 700 µL. If your volume is higher repeat this step as many times as you needed.

8. Centrifuge at 10,000 g for 30 seconds. Remove the flow-through and place back High-Q™ Spin Column into a Collection Tube.
9. Add 700 µL BQ2 Buffer.
 - ✓ Check absolute ethanol has been added to BQ2 Buffer.

10. Centrifuge at 10,000 g for 30 seconds. Remove the flow-through and place back the High-Q™ Spin Column into a Collection Tube.

11. To dry High-Q™ Spin Column, centrifuge at 10,000 g for 1 minute.

12. Place the High-Q™ Spin Column into a clean 1.5 mL Tube.


13. Add 30 - 60 µL prewarmed Elution Buffer or Water (Molecular Biology Grade) to elute purified DNA.

Prewarm Elution Buffer or Water at 70°C.

14. Incubate at room temperature, 2 minutes.

15. Centrifuge at 10,000 g for 1 minute.

16. Store DNA at -20°C.



✓ Ethanol has been added

II. DNA PURIFICATION FROM ENZYMATIC REACTION

1. Add 5 volumes of BCL1 Buffer to 1 volume of DNA solution (PCR, enzymatic reaction, etc).
2. Mix well by pipette.
3. To guarantee optimal DNA adsorption, add 3-5 μL 3M Sodium Acetate pH=5.0 and mix well.

4. Transfer the mix to a High-Q™ Spin Column placed into a Collection Tube.

The maximum loading volume of the High-Q™ Spin Column is 700 μL . If your volume is higher repeat step 3 as many times as you needed.

5. Centrifuge at 10,000 g for 30 seconds. Remove the flow-through and place back the High-Q™ Spin column into a Collection Tube.

6. Add 700 μL BQ2 Buffer.

✓ Check absolute ethanol has been added to BQ2 Buffer.

7. Centrifuge at 10,000 g for 30 seconds. Remove the flow-through and place back the Spin column into a Collection Tube.

8. To dry High-Q™ Spin Column and eliminate residual ethanol, centrifuge again at 10,000 g for 1 minute.

9. Place the High-Q™ Spin Column into a clean 1.5 mL Tube.

10. Add 30 - 60 μL prewarmed Elution Buffer or Water (Molecular Biology Grade) to elute purified DNA.

Prewarm Elution Buffer or Water at 70°C.

11. Incubate at room temperature, 2 minutes.

12. Centrifuge at 10,000 g for 1 minute.

13. Check DNA quantity by spectrophotometry and quality on agarose electrophoresis gel using COBAL™ Loading Buffer 6x provided.

14. Store at -20°C.

