High-Q™ Spin-Column DNA Cleanup Kit

Ordering info

TBK0195, 3 reactions (sample)

TBK0196, 50 reactions

TBK0197, 200 reactions

TBK0198, 250 reactions

Description

High-Q™ Spin-Column DNA Cleanup Kit is an easy silica-membrane-based DNA purification kit to remove salts, proteins, primers or other contaminants in a DNA solution. Purification of DNA embedded in reaction mixtures like PCR or any other enzymatic reaction.

Features

- Complete removal of DNA contaminants.
- Purification of small fragments ≥ 75 bp.
- High DNA recovery, $>80\%$.
- No phenol extraction.
- Easy and Fast protocol.

Applications

- DNA isolation from enzymatic reaction mixtures such as PCR, restriction digestion, labelling reactions, dephosphorylation, etc.
- DNA obtained is suitable for downstream molecular biology applications such as cloning, PCR, sequencing, digestion, genotyping, etc.

Quality Control

Isolation of a 0.5 kb DNA fragment from PCR reaction is checked by: integrity (agarose gel electrophoresis), quantity and quality ($A_{260}/A_{280} = 1.8 \pm 0.2$).

Kit Components

Components	TBK0196	TBK0197	TBK0198
High-Q™ DNA Spin Columns with Collection Tubes	50	200	250
BCL1 Buffer	30 mL	120 mL	2 x 65 mL
BQ2 Buffer	8 mL ^a	35 mL ^b	40 mL ^c
Sodium Acetate 3M, pH=5	1 mL	1 mL	1 mL
Elution Buffer	15 mL	15 mL	25 mL
COBAL™ DNA Loading Buffer 6x	0.1 mL	0.5 mL	0.5 mL

Order Info Kit Components: High-Q™ DNA Spin Columns (TBM0017) | Collection Tubes (TBM0020) | BCL1 Buffer (TBB0535) | BQ2 Buffer (TBB0529) | Elution Buffer (TBB0510) | COBAL™ DNA Loading Buffer 6x (TBB0321).

Before its use:

- ^a Add 32 mL absolute ethanol and mix well.
- ^b Add 140 mL absolute ethanol and mix well.
- ^c Add 160 mL absolute ethanol and mix well.

Storage

Store the kit at 25°C.

Store COBAL™ DNA Loading Buffer 6x at -20°C

Material required (not supplied)

- Ethanol (CAS 64-17-5).
- 1.5 mL Tubes.

PROTOCOL

1. Add 5 volumes of BCL1 Buffer to 1 volume of DNA solution (PCR, enzymatic reaction, etc).
2. Mix well by pipette.
3. [Optional] Because optimal DNA adsorption requires a $\text{pH} \leq 7.5$, in some cases it is necessary to add 3-5 μL 3M Sodium Acetate $\text{pH}=5.0$ to adjust pH.
4. Transfer the mix to a High-Q™ DNA Spin Column placed into a Collection Tube.

The maximum loading volume of the Spin Column is 700 μL . If your volume is higher repeat steps 4 and 5 as many times as you needed.


5. Centrifuge at 10,000 g for 30 seconds. Remove the flow-through and place back the High-Q™ Spin column into a Collection Tube.
6. Add 700 μL BQ2 Buffer.

✓ Check absolute ethanol has been added to BQ2 Buffer.

7. Centrifuge at 10,000 g for 30 seconds. Remove the flow-through and place back the High-Q™ Spin column into a Collection Tube.
8. To dry High-Q™ Spin Column and eliminate residual ethanol, centrifuge again at 10,000 g for 1 minute.
9. Place the High-Q™ Spin Column into a clean 1.5 mL Tube.
10. Add 30-60 μL prewarmed Elution Buffer or Water (Molecular Biology Grade) to elute purified DNA.

Prewarm Elution Buffer or Water at 70°C.

11. Incubate at room temperature, 2 minutes.
12. Centrifuge at 10,000 g for 1 minute.
13. Check DNA quantity by spectrophotometry and quality on agarose electrophoresis gel using COBALT™ Loading Buffer 6x provided.
14. Store at -20°C.

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