

# High-Q™ Spin Column TIARIZOL™ PLUS RNA Purification Kit

## Ordering info

TBK0243, 3 reactions (sample)

TBK0244, 50 reactions

TBK0245, 200 reactions

## Description

High-Q™ Spin Column TIARIZOL™ PLUS RNA Purification Kit is an ideal kit to obtain high quality RNA from a wide variety of samples, including cultured cells, animal and plant tissues, blood, bacteria and yeast. The kit combines the effective lysis of TIARIZOL™ Reagent with the efficient purification system based on silica High-Q™ RNA spin columns to provide ultrapure total RNA. The procedure includes sequential steps: lysis, followed by silica adsorption and washing steps. It is not necessary a phase separation, or post-purification steps. The kit allows to make high amounts of RNA which are very clean due to the column step.

## Features

- High yield and purity (A260/A280 ~1.8).
- RNA suitable for downstream applications.

## Applications

- cDNA library
- Northern blotting
- Nuclease protection assay
- RT-PCR
- *In vitro* translation

## Quality Control

Functionally tested with cultured cells.

## Kit Components

Components	TBK0244	TBK0245
TIARIZOL™ Reagent	50 mL	3 x 60 mL
High-Q™ RNA Spin Column	50 units	200 units
DNase I (5 U/μL)	250 μL	2x 500 μL
10x DNase-I Buffer	1.5 mL	2 x 1.5 mL
WRNA-1 Buffer	20 mL <sup>a</sup>	70 mL <sup>b</sup>
WRNA-2 Buffer	12 mL <sup>c</sup>	45 mL <sup>d</sup>
Water, nuclease free	15 mL	25 mL

**Order Info Kit Components:** TIARIZOL™ Reagent (TBR0102) | High-Q™ RNA Spin Column with Collection Tubes (TBM0012) | WRNA-1 Buffer (TBB0544) | WRNA-2 Buffer (TBB0545) | Water, nuclease free (TBB0302).

## Before its use

- <sup>a</sup> Add 12 mL of Absolute Ethanol
- <sup>b</sup> Add 42mL of Absolute Ethanol
- <sup>c</sup> Add 48 mL of Absolute Ethanol
- <sup>d</sup> Add 180 mL of Absolute Ethanol

## Storage

Store the kit at 25°C.

Store DNase-I and 10x DNase-I Buffer at -20° C.

Store TIARIZOL™ Reagent at 4°C. Protect from light.

## Reagents Required (not supplied)

- Ethanol (CAS 64-17-5).
- 1.5 mL microcentrifuge tubes
- 2 mL tubes

Also available:

TIARIZOL™ Reagent (TBR0100)

## PROTOCOL

### I. SAMPLE PREPARATION

**Cultured cells:** Harvest cells by centrifugation and decant the supernatant. Resuspend the pellet with the supernatant remaining. Add 0.8 mL of TIARIZOL™ Reagent per  $5-10 \times 10^6$  cells. Lyse cells by repetitive pipetting up and down.

**Bacterial Cells:** From a stationary culture harvest bacterial cells by centrifugation. Resuspend the pellet with 200  $\mu$ L Tris-HCl 10 mM, pH=7.5. Add 10  $\mu$ L Lysozyme (50 mg/mL) and mix by vortex. Incubate at 37°C for 30 minutes. Then add 0.8 mL of TIARIZOL™ Reagent and homogenize well.

**Yeast:** From a stationary culture harvest yeast cells by centrifugation. Disruption may require a homogenizer. Then add 0.8 mL of TIARIZOL™ Reagent and homogenize well.

**Plant Tissue:** Use up to 50 mg of fresh vegetable material and up to 20 mg of dried material. Plant material can be homogenized with liquid nitrogen using a mortar and a pestle or using any commercially available equipment for homogenization. Add 0.8 mL of TIARIZOL™ Reagent to homogenized material and mix well.

**Animal tissue:** Use up to 5 mg of animal tissue. Homogenization can be done with liquid nitrogen using a mortar and a pestle or using any commercially available equipment for homogenization. Add 0.8 mL of TIARIZOL™ Reagent to homogenized material. Mix well.

*After addition of TIARIZOL™ Reagent, incubate 5 minutes at room temperature. Samples with TIARIZOL™ Reagent can be stored at -80°C until RNA extraction procedure.*

### II. RNA EXTRACTION

1. If the samples were stored at -80°C, thaw and homogenize them.  
*A homogeneous sample is essential before RNA Purification steps.*
2. Centrifuge at 13,000 g for 5 minutes at 4°C and transfer the supernatant to a 2 mL clean tube.
3. Add **1 volume of absolute ethanol** and mix.
4. Transfer **700  $\mu$ L of the mixture** in the High-Q™ RNA Spin Column reservoir coupled to a Collection Tube.
5. Centrifuge at 10,000g for 1 minute. Discard the flow-through.
6. Repeat steps 4-5 for the remaining mixture obtained in step 3.
7. To dry silica matrix, centrifuge at 10,000g for 1 minute.
8. Add **50  $\mu$ L DNase Mixture** (5  $\mu$ L DNase-I + 45  $\mu$ L 10x DNase-I Buffer) in the center of High-Q™ RNA Spin Column.
9. Incubate for 15 minutes at room temperature (15-25°C).
10. Add **500  $\mu$ L WRNA-1 Buffer** (✓) and centrifuge at 10,000 g for 1 minute. Discard the flow-through and place the High-Q™ RNA Spin Column back into the Collection Tube.  
*✓ Check Ethanol has been added.*
11. Add **700  $\mu$ L WRNA-2 Buffer** (✓) and centrifuge at 10,000 g for 1 minute. Discard the flow-through.  
*✓ Check Ethanol has been added.*



12. Place High-Q™ RNA Spin Column back in the Collection Tube and repeat step 11.
13. To dry silica matrix, centrifuge at 13,000g for 1 minute.
14. Place High-Q™ RNA spin column into a clean 1.5 mL microcentrifuge tube (RNase-free).
15. Carefully and without touching the matrix, add in the center of High-Q™ RNA Spin Column, **50-100 µL of molecular biology grade water** (nuclease free) prewarmed at 65°C.
16. Incubate at room temperature for 1-2 minutes.
17. Centrifuge for 1 minute at 13,000 g. RNA isolated is in the eluate. Discard High-Q™ RNA Spin Column.
18. Store at -80°C.