# High-Q™ -Spin-Column Plant RNA Purification Kit

## Ordering info

TBK0278, 3 reactions (sample)

TBK0279, 20 reactions

TBK0280, 50 reactions

TBK0281, 100 reactions

## Description

High-Q™-Spin-Column Plant RNA Purification Kit is an easy silica-membrane-based system for RNA purification from a wide variety of plant species. An optimized lysis buffer guarantees a good yield while the use of High-Q™ RNA Spin Columns allows a good quality RNA, suitable for downstream applications.

## Features

- **Safety**, no phenol extraction, no ethanol precipitation.
- **High yield and purity** (A260/A280 ~2.0; A260/A230 ~2.0-2.2).
- **Isolated RNA is ready to use for downstream molecular biology applications.**
- **Easy and fast protocol.**

## Applications

- Purification of RNA from plant tissue, including plant cells, leaves, seeds, fruits, or roots.
- Purification of RNA plant using different starting plant materials: frozen, fresh, or dried.
- RNA obtained is suitable for downstream molecular biology applications such as RT-PCR, RT-qPCR, Northern, cDNA library, nuclease protection assay, *in vitro* translation, etc.

## Quality Control

RNA purified is checked by integrity (agarose gel electrophoresis), quantity and quality (A260/280 ~2.0).

## Kit Components

Components	TBK0279	TBK0280	TBK0281
High-Q™ RNA Spin Column with Collection Tubes	20	50	100
BPRNA-1 Buffer	15 mL <sup>a</sup>	35 mL <sup>a</sup>	60 mL <sup>a</sup>
BPRNA-2 Buffer	1.8 mL	2 x 1.8 mL	10 mL
DNase I (5 U/μL)	100 μL	250 μL	500 μL
10x DNase-I Buffer	1.5 mL	2 x 1.5 mL	10 mL
WRNA-1 Buffer	10 mL <sup>b</sup>	20 mL <sup>c</sup>	35 mL <sup>d</sup>
WRNA-2 Buffer	6 mL <sup>e</sup>	12 mL <sup>f</sup>	25 mL <sup>g</sup>
Water, nuclease free	5 mL	5 mL	10 mL

**Order Info Kit Components:** High-Q™ RNA Spin Columns (TBM0012) | BPRNA-1 Buffer (TBB0555) | BPRNA-2 Buffer (TBB0556) | DNase-I (TBZ0320) | 10x DNase Buffer (TBB0319) | WRNA-1 Buffer (TBB0544) | WRNA-2 Buffer (TBB0545) | Water, nuclease free (TBB0302).

**Components for samples are ready to use!**

### Before its use:

- <sup>a</sup> Add 10 μL β-mercaptoethanol per 1 mL BPRNA-1 Buffer.
- <sup>b</sup> Add 6 mL absolute ethanol and mix well.
- <sup>c</sup> Add 12 mL absolute ethanol and mix well.
- <sup>d</sup> Add 21 mL absolute ethanol and mix well.
- <sup>e</sup> Add 24 mL absolute ethanol and mix well.
- <sup>f</sup> Add 48 mL absolute ethanol and mix well.
- <sup>g</sup> Add 100 mL absolute ethanol and mix well.

## Storage


Store the kit at 25°C and DNase-I at -20°C.

## Material required (not supplied)

- 1.5 mL Microcentrifuge tubes (RNase free).
- Ethanol (CAS 64-17-5).
- β-mercaptoethanol (βME) (CAS 60-24-2)

## PROTOCOL

- Grind up to 100 mg of plant material in liquid nitrogen using a mortar and a pestle. With a freeze spatula, collect the powder into a frozen 1.5 mL tube.
- Add **600  $\mu$ L BPRNA-1 Buffer** and mix vigorously by vortex.  
*Check  $\beta$ ME has been added. All vegetal powder must be mixed with the buffer.*
- Add **60  $\mu$ L BPRNA-2 Buffer** and mix vigorously by vortex.
- Incubate at 50 °C for 5 minutes. Mix by inversion from time to time.
- Centrifuge at 13,000 g for 5 minutes.
- With a pipette, transfer the supernatant very carefully to another tube, avoiding vegetal debris.
- Add **0.5 volumes of absolute ethanol** (~ 250  $\mu$ L). Mix by inversion.
- Transfer up 700  $\mu$ L mixture to a High-Q™ RNA Spin Column placed into a Collection Tube.
- Centrifuge at 10,000 g for 1 minute. Remove the flow-through and place back the High-Q™ RNA Spin Column into a Collection Tube. If necessary, repeat steps 8 and 9 with the remaining mixture.
- Centrifuge at 10,000 g for 1 minute to dry the column matrix.
- Add **50  $\mu$ L DNase Mixture** in the center of High-Q™ RNA Spin Column.  
*DNase Mixture: Mix with pipette 5  $\mu$ L DNase-I + 45  $\mu$ L 10x DNase-I Buffer. Avoid vortex.*
- Incubate for 15 minutes at room temperature (15-25°C).
- Add **500  $\mu$ L WRNA-1 Buffer** (✓) and centrifuge at 10,000 g for 1 minute. Discard the flow-through and place the High-Q™ RNA Spin Column back into the Collection Tube.  
✓ Check Ethanol has been added.
- Add **500  $\mu$ L WRNA-2 Buffer** (✓).  
✓ Check Ethanol has been added.
- Centrifuge at 10,000 g for 1 minute. Discard flow-through. Place High-Q™ RNA Spin Column back in the Collection Tube and repeat from step 14.
- To dry silica matrix, centrifuge at 10,000g for 1 minute.
- Place High-Q™ RNA spin column into a clean 1.5-mL microcentrifuge tube.
- Carefully and without touching the matrix, add in the center of High-Q™ RNA Spin Column, **50-100  $\mu$ L Water** (nuclease free) prewarmed at 65°C.
- Incubate at room temperature for 1-2 minutes.
- Centrifuge for 1 minute at 13,000 g. RNA isolated is in the eluate. Discard High-Q™ RNA Spin Column.
- Store at -80°C.

# High-Q™-Spin-Column Plant RNA Purification Kit

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