Tiaris™- One-Step RT-PCR Kit

Ordering Info

TBK1040, 100 rxns of 25 µL

Description

Tiaris™- One-Step RT-PCR Kit is designed for first-strand cDNA synthesis and subsequent PCR in a single-tube reaction procedure, decreasing contamination risk and reducing hands-on time considerably. The kit consists of a RT mix and a PCR Master Mix. The RT mix comprises a modified M-MuLV reverse transcriptase with increased thermostability and lacking RNase H activity, and an advanced RNase inhibitor to impede RNA degradation. The Master mix contains all the reagent (except PCR primers and template) needed for running PCR reactions.

Tiaris™- One-Step RT-PCR Kit can be used with any RNA template, including mRNA, viral RNA and total RNA, as the RTase is not inhibited by rRNA or tRNA. The kit is optimized for rapid one-step, one-tube RNA quantification, reducing experimental variation and contamination with a convenient RT-PCR protocol.

Features

- Thermostable Reverse Transcriptase, allows for cDNA synthesis at 45°C-55°C.
- RNase H deficient.
- Fast Hot-Start PCR Master Mix.
- High yield, specificity, and sensitivity.

Kit Components

Components	TBK1040
PCR Master Mix 2x*	1.25 mL
RT Mix	125 µL
PCR Grade Water, nuclease free	2 mL

* The PCR Master Mix contains the hot-start Taq Polymerase, 6mM MgCl₂, 2mM dNTPs, enhancers and stabilizers.

Storage

Shipped on blue ice. Upon receipt, kit components should be immediately stored at -20°C. Avoid repeated freezing and thawing. Maintain cold when thawed.

Applications

- One-step RT-PCR.
- Gene expression analysis by endpoint RT-PCR.

Technical Assistance

Please refer any technical questions to support@tiarisbiosciences.com

PROTOCOL

1. Gently vortex and briefly centrifuge kit components after thawing.
2. Place a tube on ice and add the following components for each 20 μ L reaction. Prepare sufficient master mix for the number of reactions. Consider one or two extras:

Components	Volume
PCR Master Mix (2x)	12.5 μ L
Forward primer (10 μ M)	1 μ L
Reverse primer (10 μ M)	1 μ L
RT Mix	1.25 μ L
Template RNA (step 5)	*
PCR Grade Water, nuclease-free	up 25 μ L
Final Volume	25 μL

3. Dispense the master mix into wells of PCR plate.
4. Gently vortex and spin down the samples.
5. Add in each well the RNA sample. Mix well by pipetting.
Template RNA: in case of total RNA 1 pg-1 μ g, in case of mRNA > 0.01 pg (the RNA sample should be completely free of contaminating genomic DNA).
6. Seal the PCR plate with optical film.
7. Set-up RT-PCR cycling:

Suggested thermal cycling conditions:

Process	Cycles	Temperature	Time
cDNA synthesis	1 x	45°C *	10-15 min
Enzyme Activation	1 x	95 °C	3 min
Denaturation	30-40 x	95 °C	10 sec
Annealing		Ta **	10 sec
Extension	1 x	72 °C	15 sec/ kb
Final Extension	1 x	72 °C	1 min
Hold	1 x	4 °C	∞

* For most applications, cDNA synthesis should be carried out at 45°C. In case of targets with high secondary structure, synthesis may be optimized by carrying out at 55°C.

**Ta depends on the Tm of both primers, which in turn depend on the buffer composition. Ideally, primers should have a predicted Tm of about 60°C. In general, using a Ta of 1-3°C above the Tm of the primer with the lowest Tm is a good starting point. It is highly recommended to determine optimum annealing temperature experimentally by performing a temperature gradient.

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