

# LGex 3410

## MediLine



<b>Total gross capacity</b>	310 litre
<b>Ext. dimensions</b> in mm (w/d/h)	600/615/1840
<b>Int. dimensions</b> in mm (w/d/h)	420/400/1587
<b>Energy consumption</b> in 24 h	1,309 kWh
Klimaklasse <sup>1</sup>	SN-ST
Cooling system	static
Defrosting	manual
Temperature range	-9°C to -30°C
Gastronorm	-
Housing material/colour	steel/white
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	8
Adjustable Shelving	-
Number of baskets	-
Maximum number of baskets	8
Shelf material	evaporator plates
Shelf loading	24 kg
Height of adjustable legs	-
Castors/type of castors	-
Door handle	Ergonomic door handle
Foot pedal opener	-
Self-closing door	yes
Insulation	7 - 7 mm
Defrost water drain	-
Refrigerant	R 600a
Noise level	-
Gross/net weight	94/88 kg
Interior light	-
EAN-Nr.	9005382136515

Climate Ratings: **SN** +10°C to +32°C    **N** +16°C to +32°C    **ST** +16°C to +38°C    **T** +16°C to +43°C